GET STARTED

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IT I PRIMI PASSI CON KOLIBREE

KLTB001

IMPORTANT SAFEGUARDS

This appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved. Children should not play with the appliance. Cleaning and user maintenance should not be made by children without supervision. The appliance is only to be used with the power supply unit provided with the appliance. The external flexible cable or cord of this adaptor cannot be replaced; if the cord is damaged, the adaptor shall be scrapped. The socket-outlet should be installed near the equipment and be easily accessible.

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


ES ENCENDER EL CEPILLO DE DIENTES
IT ACCENDERE LO SPAZZOLINO DA DENTI




Press the on/off button, the LED blinks briefly before turning on.

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Press and hold the button until the blue LED starts to blink.

FR Maintenez le bouton appuyé jusqu'à ce que la lumière clignote en bleu.
ES Mantenga pulsado el botón hasta que la luz parpadee en azul.
DE Halten Sie die Taste gedrückt, bis die Kontrolllampe blau blinkt.
IT Tenere premuto il pulsante fino a che la luce blu lampeggia.

The toothbrush is made discoverable to Bluetooth® for 60 seconds.

FR La brosse est détectable par Bluetooth® pendant 60 secondes.
ES El cepillo se puede detectar por Bluetooth® durante 60 segundos.
DE Die Zahnbürste ist 60 Sekunden lang per Bluetooth® erkennbar.
IT Lo spazzolino è individuabile tramite Bluetooth® per 60 secondi.

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

Download and install the Kolibree application on the App Store.

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Download and install the Kolibree application on the Play Store.

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

Create your profile on the application, then choose a game option to play while brushing your teeth.

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









Turn the used head until the two indicators are aligned.

FR Tournez la brosse jusqu'à ce que les deux indicateurs soient alignés.
ES Gire el cabezal hasta que se alineen los dos indicadores.
DE Drehen Sie den Bürstenkopf bis die beiden Markierungen einander gegenüberliegen.
IT Ruotare la testina fino a che i due indicatori sono allineati.

Pull the brush head up and replace it with a new one.

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WWW.KOLIBREE.COM

Order your brush heads at any time on the Kolibree website!

FR Commandez vos brosses à tout moment sur le site Kolibree !
ES ¡Adquiera sus cabezales en cualquier momento en la página de Kolibree!
DE Neue Bürstenköpfe können Sie jederzeit auf der Kolibree-Website bestellen!
IT Ordinare le testine nuove nel sito Kolibree!

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
Simply place the toothbrush on the base for charging.

FR Posez simplement la brosse à dents sur son socle pour la recharger.
ES Simplemente coloque el cepillo de dientes en su base para recargarlo.
DE Um die Zahnbürste aufzuladen, stellen Sie sie einfach auf den Ladesockel.
IT Per ricaricarla è sufficiente appoggiare lo spazzolino da denti sulla base.

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SWITCHING MODE POWER SUPPLY


PLEASE READ THESE INSTRUCTIONS CAREFULLY BEFORE USING THIS PRODUCT :

The external flexible power cord of this transformer cannot be replaced; if the cord is damaged in any way, the transformer should be replaced.
In the event of a short circuit or power overload, the protective devices cannot be reset or replaced, and the product should be replaced.
Electrical function of the transformer : Non-inherently short circuit proof safety isolating transformer.
Kolibree's transformer is only intended for indoor use.
If the plug on the supply cord is serving as a disconnection device, the socket outlet shall be installed near the Kolibree base and shall be easily accessible.
If the mains plug is used as the disconnection device, the disconnection device shall remain readily operable.
The battery must be removed from the appliance before it is scrapped.
The appliance must be disconnected from the supply mains when removing the battery.
The battery is to be disposed of safely.

CAUTION: RISK OF ELECTRIC SHOCK
Dry location use only.
Disconnect power to the receptacle before installing or removing the unit.
When removing receptacle cover screws, cover may fall across plug pins or receptacle may become displaced.
Use only with duplex receptacle having center screw.
Secure unit in place by receptacle cover screw.


REMOVING BATTERY

Please note this process is not reversible.

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ENGLISH

IMPORTANT PRODUCT INFORMATION

For safety, compliance, recycling, and other important information regarding your Kolibree toothbrush, please refer to the Kolibree toothbrush user's guide which is accessible on your Kolibree application.

KOLIBREE: LICENSE AGREEMENT AND TERMS OF USE

BEFORE USING THE KOLIBREE TOOTHBRUSH DEVICE, PLEASE READ THE LICENSE AGREEMENT AND TERMS OF USE FOUND AT WWW.KOLIBREE.COM AND ALL RULES AND POLICIES FOR THE KOLIBREE TOOTHBRUSH DEVICE AND SERVICES RELATED TO THE DEVICE (INCLUDING, BUT NOT LIMITED TO, THE KOLIBREE PRIVACY NOTICE AND ANY RULES OR USAGE PROVISIONS ACCESSIBLE THROUGH THE WWW.KOLIBREE.COM WEBSITE OR THE KOLIBREE STORE) (COLLECTIVELY, THIS "AGREEMENT"). BY USING THE KOLIBREE TOOTHBRUSH DEVICE YOU AGREE TO BE BOUND BY THE TERMS OF THIS AGREEMENT.

Contact Information. For communications concerning this Agreement, you may contact Kolibree by writing to Kolibree SAS, 36 rue de Turin, 75008 Paris, France.

LIMITED WARRANTY

This Warranty for the Kolibree toothbrush Device (the "Device") is provided by Kolibree

SAS, 36 rue de Turin, 75008 Paris, France.

We warrant the Device against defects in materials and workmanship under ordinary consumer use for a period of one year or the current legal period after the initial date and place of purchase. During this warranty period, if a defect arises in the Device, and you follow the instructions for returning the Device, we will, at our option and to the extent permitted by law, either (i) repair the Device using either new or refurbished parts, (ii) replace the Device with a new or refurbished Device, or (iii) refund to you all or part of the purchase price of the Device. This limited warranty applies to any repair, replacement part or replacement Device for the remainder of the original warranty period or for ninety days, whichever period is longer. All replaced parts and Devices for which a refund is given shall become our property. This limited warranty applies only to hardware components of the Device that are not subject to accident, misuse, neglect, fire or other external causes, alterations, repair, or commercial use.

Instructions. For specific instructions about how to obtain warranty service for your Device, please contact Customer Service using the contact information provided in this guide. In general, you will need to deliver your Device in either its original packaging or in equally protective packaging to the address specified by Customer Service. Before you deliver your Device for warranty service, it is your responsibility to back up any data, software, or other materials you may have stored or preserved on your Device. It is possible that such data, software or other materials will

be lost or reformatted during service, and we will not be responsible for any such damage or loss.

Limitations. TO THE EXTENT PERMITTED BY LAW, THE WARRANTY AND REMEDIES SET FORTH ABOVE ARE EXCLUSIVE AND IN LIEU OF ALL OTHER WARRANTIES AND REMEDIES, AND WE SPECIFICALLY DISCLAIM ALL STATUTORY OR IMPLIED WARRANTIES, INCLUDING, BUT NOT LIMITED TO, WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND AGAINST HIDDEN OR LATENT DEFECTS. IF WE CANNOT LAWFULLY DISCLAIM STATUTORY OR IMPLIED WARRANTIES, THEN TO THE EXTENT PERMITTED BY LAW, ALL SUCH WARRANTIES SHALL BE LIMITED IN DURATION TO THE DURATION OF THIS EXPRESS LIMITED WARRANTY AND TO REPAIR OR REPLACEMENT SERVICE.

SOME JURISDICTIONS DO NOT ALLOW LIMITATIONS ON HOW LONG A STATUTORY OR IMPLIED WARRANTY LASTS, SO THE ABOVE LIMITATION MAY NOT APPLY TO YOU. WE ARE NOT RESPONSIBLE FOR DIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM ANY BREACH OF WARRANTY OR UNDER ANY OTHER LEGAL THEORY. IN SOME JURISDICTIONS THE FOREGOING LIMITATION DOES NOT APPLY TO DEATH OR PERSONAL INJURY CLAIMS, OR ANY STATUTORY LIABILITY FOR INTENTIONAL AND GROSS NEGLIGENT ACTS AND/OR OMISSIONS, SO THE ABOVE EXCLUSION OR LIMITATION MAY NOT APPLY TO YOU. SOME JURISDICTIONS DO NOT ALLOW THE EXCLUSION OR LIMITATION OF

DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES, SO THE ABOVE EXCLUSION OR LIMITATION MAY NOT APPLY TO YOU. THIS "LIMITATIONS" SECTION DOES NOT APPLY TO CUSTOMERS IN THE UNITED KINGDOM AND THE EUROPEAN UNION.

This limited warranty gives you specific legal rights. You may have additional rights which vary from jurisdiction to jurisdiction, and this limited warranty does not affect those rights.

You can see a version of the Kolibree License Agreement and Terms of Use, One-Year Limited Warranty and other Kolibree information in other languages at <http://www.kolibree.com>.

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FCC STATEMENT

Warning: Changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications.

However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures: Reorient or relocate the receiving antenna. Increase the separation between the equipment and receiver. Connect the equipment into an outlet on a circuit different from that to which the receiver is connected. Consult the dealer or an experienced radio/TV technician for help.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

IC STATEMENT

This device complies with Industry Canada licence-exempt RSS standard(s). Operation is subject to the following two conditions: (1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device.

Under Industry Canada regulations, this radio transmitter may only operate using an antenna of a type and maximum (or lesser) gain approved for the transmitter by Industry Canada. To reduce potential radio interference to other users, the antenna type and its gain should be so chosen that the equivalent isotropically radiated power (e.i.r.p.) is not more than that necessary for successful communication.

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Designed by Kolibree in France. Assembled in the Philippines.

Kolibree © 2014 All rights reserved. iPhone, iPod and iPad are registered trademarks of Apple Inc, registered in the US and other countries. Android is a trademark of Google Inc.

 @Kolibree

WWW.KOLIBREE.COM

 /sharekolibree